


i3 GOGGLES



The all-new *i3* goggles* are offered in 2 models...
monocular and binocular.

Both versions come with a pre-attached hat clip,
designed to be clipped to the bill of a hat or glasses.

A head strap is also included with the binocular
version, along with hardware for mounting.

All models are include:

- 2 sets of rubber glare shields, one with a flare to shield sunlight and one without for better peripheral vision
- user-replaceable cables

Both models are small and lightweight, and are
compatible with Ibex EVO, Pro, and Lite as well as
with Bantam and Easi Scan 1-4.



Custom designed optics

- 800x600 color Oled display(s)
- half the weight of InSite2 goggles
- specially designed eye cups for inside and outdoor viewing
- adjustable links, universal clip
- GoPro mount
- 3 year warranty [regular use]



**patent pending*

For more information or to schedule a demonstration,
contact us at 866.365.6596 • info@eimedical.com

E.I. Medical Imaging® | *ibex*®